



ARROW



ARROW-TREROTOLA™ PTD®

Fast and easy removal of thrombi



ORDERING INFORMATION

ARROW-TREROTOLA PTD 5 FR.						ARROW
REF.	CATHETER LENGTH	FRAGMENTATION BASKET	SHEATH INCLUDED	ROTATOR DRIVE UNIT	QTY	
PT-03000-R	–	–	–	3000 rpm	1/case	
PT-65509*	65 cm	9 mm	–	–	1/case	
PT-45509**	65 cm	9 mm	2/6 Fr.	–	1/case	
PT-65509-HFC	65 cm	9 mm	2/6 (HF) Fr.	3000 rpm	1/case	

ARROW-TREROTOLA OVER-THE-WIRE PTD 7 FR.						ARROW
REF.	CATHETER LENGTH	FRAGMENTATION BASKET	MAX. GUIDEWIRE COMPATIBILITY	SHEATH INCLUDED	ROTATOR DRIVE UNIT	QTY
PT-03009-RW	–	–	–	–	3000 rpm	1/case
PT-65709-W***	65 cm	9 mm	0.025"	–	–	1/case
PT-65709-WC****	65 cm	9 mm	0.025"	2/7 Fr.	3000 rpm	1/case
PT-65709-HFWC****	65 cm	9 mm	0.025"	2/7 (HF) Fr.	3000 rpm	1/case
PT-12709-WC****	120 cm	9 mm	0.025"	2/7 Fr.	3000 rpm	1/case

PTD ACCESSORY COMPONENTS INTRODUCER SHEATHS								ARROW
REF.	SHEATH SIZE	SHEATH LENGTH	VESSEL DILATOR LENGTH	MAX. GUIDEWIRE COMPATIBILITY	RADIOPAQUE TIP MARKER	LARGE-BORE-SIDEARM	COLOUR-CODED HUB	QTY
CL-08505	5 Fr.	2"	5"	0.038"	–	–	grey	10/case
CL-08605	6 Fr.	2"	5"	0.038"	–	–	green	10/case
CL-08605-HF	6 Fr.	2"	5"	0.038"	•	•	green	5/case
CL-08705-HF	7 Fr.	2"	5"	0.038"	•	•	orange	5/case

Does not contain natural rubber latex.

* When ordering this component, the PT-03000-R and CL-08605-HF must also be ordered.

** When ordering this component, the PT-03000-R must also be ordered.

*** When ordering this component, the PT-03009-RW and CL-08705-HF must also be ordered.

**** Each product includes: two radiopaque polyurethane sheaths with integral side port/hemostasis valve, and one vessel dilator with SnapLock™ feature.

LIBRARY OF CLINICAL SUPPORT

Hein, A.N., Vesely, T.M. Use of the Percutaneous Thrombolytic Device for the Treatment of Thrombosed Pseudoaneurysms During Mechanical Thrombectomy of Hemodialysis Grafts. *Journal of Vascular Interventional Radiology*. 2002; 13: 201–204.

Lajvardi, A., Trerotola, S.O., Strandberg, J.D., Samphilipo, M.A., Magee, C. Evaluation of Venous Injury Caused by a Percutaneous Mechanical Thrombolytic Device. *Cardiovascular Interventional Radiology*. 1995; 18: 172–178.

Lazzaro, C.R., Trerotola, S.O., Shah, H., Namyslowski, J., Moresco, K., Patel, N. Modified Use of the Arrow-Trerotola Percutaneous Thrombolytic Device for the Treatment of Thrombosed Hemodialysis Access Grafts. *Journal of Vascular Interventional Radiology*. 1999; 10: 1015–1031.

McLennan, G., Trerotola, S.O., Davidson, D. et al. The Effects of a Mechanical Thrombolytic Device on Normal Canine Vein Valves. *Journal of Vascular Interventional Radiology*. 2001; 12:89–94.

Patel AA, Tuite CM, Trerotola SO. *Cardiovasc. Intervent Radiol*. 2005 ;28(6): 704–13 Mechanical Thrombectomy of Hemodialysis Fistulae and Grafts

Rocek, M., Peregrin, J.H., Lasovickova, J., Krajickova, D., Slaviokova, M. Mechanical Thrombolysis of Thrombosed Hemodialysis Native Fistulas With Use of the Arrow-Trerotola Percutaneous Thrombolytic Device: Our Preliminary Experience. *Journal of Vascular Interventional Radiology*. 2000; 11: 1153–1158.

Shatsky et al. *J.Vasc.Interv.Radiol*. 2005 ; 1605–11 Single-Center Experience with the Arrow-Trerotola Percutaneous Thrombectomy Device in the Management of Thrombosed Native Dialysis Fistulas.

Trerotola, S.O., Davidson, D.D., Filo, R.S., Dreesen, R.G., Forney, M. Preclinical In Vivo Testing of a Rotational Mechanical Thrombolytic Device. *Journal of Vascular Interventional Radiology*. 1996; 7: 717–723.

Trerotola, S.O., Vesely, T.M., Lund, G.B., Soulen, M.C., Ehrman, K.O., Cardella, J.F. Treatment of Thrombosed Hemodialysis Access Grafts: Arrow-Trerotola Percutaneous Thrombolytic Device Versus Pulse-Spray Thrombolysis. *Arrow-Trerotola Percutaneous Thrombolytic Device Clinical Trial*. *Radiology*. 1998; 206: 403–414.

Trerotola, S.O., McLennan, G., Davidson D., et al. Pre-clinical In Vivo Testing of the Arrow-Trerotola Percutaneous Thrombolytic Device for Venous Thrombosis. *Journal of Vascular Interventional Radiology*. 2001; 12: 95–103.

Trerotola, S.O., McLennan, G., Eclavea, A.C. et al. Mechanical Thrombolysis of Venous Thrombosis in an Animal Model With Use of Temporary Caval Filtration. *Journal of Vascular Interventional Radiology*. 2001; 12: 1075–1085.

Trerotola, S.O., Johnson, M.S., Schauwecker, D.S. et al. Pulmonary Emboli From Pulse-Spray and Mechanical Thrombolysis: Evaluation With an Animal Dialysis-Graft Model. *Radiology*. 1996; 200: 169–176.

Vesely, T.M., Hovsepian, D.M., Darcy, M.D., Brown, D.B., Pilgram, T. K. Angioscopic Observations After Percutaneous Thrombectomy of Thrombosed Hemodialysis Grafts. *Journal of Vascular Interventional Radiology*. 2000; 11: 971–977.

Vogel, P.M., Bansal, V., Marshall, M.W. Thrombosed Hemodialysis Grafts: Lyse and Wait With Tissue Plasminogen Activator or Urokinase Compared to Mechanical Thrombolysis With the Arrow-Trerotola Percutaneous Thrombolytic Device. *Journal of Vascular Interventional Radiology*. 2001; 12: 1157–1165.

ARROW-TREROTOLA™ PTD®

PERCUTANEOUS THROMBOLYTIC DEVICE

THE MOST PROVEN MECHANICAL THROMBECTOMY TECHNOLOGY AVAILABLE FOR DIALYSIS ACCESS

Innovative dialysis access products from Teleflex are designed to meet the needs of clinicians and inspired by them. The Arrow-Trerotola PTD improves your ability to achieve fast and easy removal of thrombi from both dialysis native AV fistulae and synthetic grafts.

SIMPLE

The Arrow-Trerotola PTD is easy to assemble and operate, and as it comes with a mechanical thrombectomy catheter, a hand-held disposable rotator drive unit and an introducer sheath, there is no need to purchase additional equipment.

- catheter lumen sidearm: permits catheter flushing during preparation and use
- activated spinning basket macerates the thrombus
- introducer sheath and large-bore sidearm assembly simplifies thrombus removal
- short procedure time patient can immediately return to dialysis²

SAFE

The Arrow-Trerotola PTD conforms to AV graft and AV fistulae walls, enabling you to remove clots with minimal risk of wall damage, and without the use of thrombolytics⁵.

REFERENCES

- 1 NNKF-K/DOQI Clinical Practice Guidelines for Vascular Access: Update 2000. *American Journal of Kidney Disease*. 2000; 37: 137-181.
- 2 Trerotola, S.O., Vesely, T.M., Lund, G.B., Soulen, M.C., Ehrman, K.O., Cardella, J.F. Treatment of Thrombosed Hemodialysis Access Grafts: Arrow-Trerotola Percutaneous Thrombolytic Device Versus Pulse-Spray Thrombolysis. *Arrow-Trerotola Percutaneous Thrombolytic Device Clinical Trial*. *Radiology*. 1998; 206: 403-414.

- 3 Lazzaro, C.R., Trerotola, S.O., Shah, H., Namyslowski, J., Moresco, K., Patel, N. Modified Use of the Arrow-Trerotola Percutaneous Thrombolytic Device for the Treatment of Thrombosed Hemodialysis Access Grafts. *Journal of Vascular Interventional Radiology*. 1999; 10: 1025-1031.
- 4 Rocek, M., et. al. Mechanical Thrombolysis of Thrombosed Hemodialysis Native Fistulas With Use of the Arrow-Trerotola Percutaneous Thrombolytic Device: Our Preliminary Experience. *Journal of Vascular Interventional Radiology*. 2000; 11: 1153-1158.

PROVEN

Satisfies key benchmarks to achieve better outcomes – including K/DOQI guidelines for both immediate and three-month patency and functionality rates. Arrow-Trerotola PTD is the only mechanical thrombectomy device indicated to pull the arterial plug³.

DIALYSIS ACCESS


STUDIES	TECHNICAL PATENCY RATE	3-MONTH PRIMARY RATE
K/DOQI Guidelines ¹	85%	40%
PDT vs PSPMT ²	95%	39%
Modified technique ³	100%	42%
Use of PTD and native fistulae ⁴	100%	70%

AV SYNTHETIC GRAFT



PTD device in use in a forearm-loop synthetic graft.

- 5 Lajvardi, A., Trerotola, S.O., Strandberg, J.D., Samphilipo, M.A., Magee, C. Evaluation of Venous Injury Caused by a Percutaneous Mechanical Thrombolytic Device. *Cardiovascular Interventional Radiology*. 1995; 18: 172-178.
- 6 Trerotola SO, Johnson MS, Shah H, Namyslowski J. Backbleeding technique for treatment of arterial emboli resulting from dialysis graft thrombolysis. *JVIR*. 1998;9:141-143.



**ARROW-TREROTOLA PTD
OVER-THE-WIRE (7 FR.)**

**ARROW-TREROTOLA
PTD (5 FR.)**

CATHETER LUMEN SIDEARM

- permits catheter flushing during preparation and use

AV FISTULAE



Pre-thrombectomy assessment of dialysis fistulae.



Thrombectomy using the over-the-wire PTD via access sheath.



Post-PTD thrombectomy fistulogram showing patent fistulae.

Courtesy of Miloslav Rocek, MD, PhD, Institute for Clinical & Interventional Radiology (IKEM), Dept. of Diagnostic & Interventional Radiology Videnska 800, 14000 Prague 4 Czech Republic

FOR USE IN APPLICATIONS

Dialysis AV fistulae and synthetic grafts



ACTIVATED SPINNING BASKET

- macerates the thrombus



SOFT, FLEXIBLE TIP

- designed to manoeuvre through vessel easily



**UNIQUE EXPANDABLE
9MM FRAGMENTATION
BASKET**

- conforms to different diameter walls
- shown to easily remove residual thrombus from dialysis vessel walls⁶

Teleflex is a leading global provider of specialty medical devices used for diagnostic and therapeutic procedures in critical care, urology and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety. We specialise in devices for general and regional anaesthesia, cardiac care, respiratory care, urology, vascular access and surgery and we serve healthcare providers in more than 140 countries. Teleflex also provides specialty products for medical device manufacturers.

Our well known brands include ARROW®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, LMA™, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, TAUT®, TFX OEM®, VASONOVA™ and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated.

Teleflex global operations: Australia, Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

YOUR INTERNATIONAL CONTACTS:

TELEFLEX HEADQUARTER INTERNATIONAL, IRELAND

Teleflex Medical Europe Ltd., IDA Business and Technology Park,
Dublin Road, Athlone, Co Westmeath
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73
orders.intl@teleflex.com

AUSTRALIA/NEW ZEALAND +61 (0)3 9081 0600

AUSTRIA +43 (0)1 402 47 72

BELGIUM +32 (0)2 333 24 60

CHINA (SHANGHAI) +86 (0)21 6381 6399

CHINA (BEIJING) +86 (0)10 6418 5699

CZECH REPUBLIC +420 (0)495 759 111

FRANCE +33 (0)5 62 18 79 40

GERMANY +49 (0)7151 406 0

GREECE +30 210 67 77 717

INDIA +91 (0)44-2836 5040

ITALY +39 0362 58 911

JAPAN +81 (0)3 3379 1511

NETHERLANDS +31 (0)88 00 215 00

PORTUGAL +351 22 541 90 85

SINGAPORE +65 6447 5022

SLOVAK REPUBLIC +421 (0)3377 254 28

SOUTH AFRICA +27 (0)11 807 4887

SPAIN +34 918 300 451

SWITZERLAND +41 (0)31 818 40 90

UNITED KINGDOM +44 (0)1494 53 27 61

For detailed information see www.teleflex.com

The products in this catalogue may not be available in all countries.
Please contact your local representative. All data current at time of printing (11/2013).
Subject to technical changes without further notice.

94 11 11 - 00 00 01 - REV A · MC/WM · 11 13 01