



RÜSCH S-BAG CLOSED URINE DRAINAGE SYSTEM

SAFE AND HYGIENIC: FOR MAXIMUM PROTECTION OF THE PATIENT AGAINST NOSOCOMIAL URINARY TRACT INFECTIONS

About 90% of all nosocomial urinary tract infections are caused by urinary drainage. Apart from the health risk to the patients, such infections also present huge cost implications to the hospital. Rüsçh offers a variety of products which reduce the risk of infection.

Optimum pressure balance and excellent flow performance

Anti-bacterial and leak proof ventilation: at the drip chamber and the urine drainage bag. Prevents collection of urine in the tube and the associated risk of infection.

The closed urine drainage system represents one of the most important elements of prophylaxis against infections.

Rüsçh urine drainage systems meet all the hygienic demands which are required for a high-quality closed system for long-term medical care.



AT A GLANCE

- Needle free sampling port
- Drip chamber with integrated check valve
- Anti-bacterial and leak proof ventilation
- Clear graduation in ml
- Flexible and kink-resistant tube
- Optimum fastening of the system
- Non drip base outlet
- Latex-free



Vented drip chamber with integrated check valve

The vented drip chamber promotes optimum drip collection with non reflux valve hence preventing retrograde infection. Optimal pressure compensation and safe outlet performance thanks to anti-bacterial ventilation.
Safe even under high pressure load.



Needle free sample port – safe and easy
Urine sampling is carried out at the ridged universal connector with a Luer syringe without a needle. Needle stick injuries are thus eliminated and the risk of invading bacteria is greatly reduced.



Hygienic and quick emptying

Non drip base outlet with one-hand clamp. Doghouse holder to fix the outlet tube, to prevent contamination, and to achieve additional floor clearance.

Universal ridged connector with tapered cone-like steps – suitable for all common catheters.

Safe, centred fixation with double hook. Additional safety thanks to extra cord with loop.



PRODUCT INFORMATION

Closed urine drainage system for standard drainage under long-term conditions.

Rüsch S-Bag			
Ref.	Capacity	Length of tube	Qty
67 30 44	2000 ml	120 cm	20
Sterile, single use, latex-free			

GERMANY • Teleflex Medical GmbH • P.O. Box 11 80 • 71385 Kernen
Phone +49 (0)7151 406 0 • Fax +49 (0)7151 406 150 • info.de@teleflexmedical.com
UNITED KINGDOM • Teleflex Medical • High Wycombe • Bucks. HP12 3ST
Phone +44 (0)1494 53 27 61 • Fax +44 (0)1494 52 46 50 • info.uk@teleflexmedical.com

www.ruesch.de • www.teleflex.com