



ARROW 3 LUMEN & 4 LUMEN 11 CM CENTRAL VENOUS CATHETERS

EACH CV SET CONTAINS:

central venous catheter radiopaque polyurethane with Blue FlexTip, tissue dilator, guidewire on the ARROW advancer, introducer needle 18 Ga. x 2 1/2" (6.35 cm), 5 ml luer slip syringe and catheter fixation.

EACH BR SET ADDS:

four high flow stopcocks, catheter over needle assembly 18 Ga. x 2 1/2" (6.35 cm) over 20 Ga. RW introducer needle and No 11 scalpel.

CV SETS					ARROW
REF.	LUMEN	SIZE	MATERIAL	QTY	
CV-10753*	3	7.5 Fr 11 cm	PUR	10	
CV-20753*	3	7.5 Fr 11 cm	PUR A+g	10	
CV-10854	4	8.5 Fr 11 cm	PUR	10	
CV-20854	4	8.5 Fr 11 cm	PUR A+g	10	

BR SETS					ARROW
REF.	LUMEN	SIZE	COATING	QTY	
BR-10753*	3	7.5 Fr 11 cm	PUR	10	
BR-20753*	3	7.5 Fr 11 cm	PUR A+g	10	
BR-10854	4	8.5 Fr 11 cm	PUR	10	
BR-20854	4	8.5 Fr 11 cm	PUR A+g	10	

A+g = ARROWgard Blue® CVC antimicrobial impregnation with chlorhexidine/silver sulfadiazine

PUR = Polyurethane

* introducer needle with echogenic tip

FLOW RATES AND PRIMING VOLUMES

LUMEN DESIGN	LUMEN	FLOW RATE (ML/HR)		PRIMING VOLUME (ML)	
		PUR	PUR A+g	PUR	A+g
	Distal 16 Ga	3.900	3.770	0.37	0.36
	Medial 14 Ga	7.200	6.350	0.68	0.67
	Proximal 16 Ga	2.900	2.650	0.38	0.37
	Distal 16 Ga	3.690	3.500	0.37	0.35
	Medial 1 14 Ga	6.560	8.400	0.43	0.39
	Medial 2 18 Ga	1.010	2.900	0.35	0.33
	Proximal 18 Ga	1.140	2.800	0.38	0.37



DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND · IDA Business Park · Athlone · Co Westmeath

Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflexmedical.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflexmedical.com

WWW.TELEFLEXMEDICAL.COM

94 10 86-00 00 01 · 08 10 01 · REV A · All data current at time of printing. Subject to technical changes without further notice.