

MDH

Zhuhai Mindhao Medical Technology Co., Ltd.

SLIM IS BIG



Video Bronchoscope
NEW A41






High Configuration



Multiple Size



One monitor is compatible with several operation parts



Multiple Disinfection Methods

Disinfection Method	YES/NO
ETO Gas	YES
Cidex OPA、 Glutaraldehyde	YES
Low Temperature Plasma (STERIS/STERRAD)	YES
Autoclave	NO
High Temperature Plasma	NO



MDH SLIM IS BIG

Zhuhai Mindhao Medical Technology Co., Ltd.

Address: 4/F, Phase 2, No.2 Chuangxin Si Road, Tangjiawan Town, Hi-tech Zone, Zhuhai 519085, P.R.China

TEL : 0756-3635312 FAX : 0756-3635312 Email:international@mindhao.com Website:www.mindhao.com